

ErgoStyle Flexion Distraction With Auto-Cocking - ES5823

Standard Features

- Electric Elevation from 19" to 32" (48 cm - 81 cm) with single foot control and expanding covers over lifting mechanism
- Speed: 19" to 32" (48 cm - 81 cm) in 13 seconds
- Manual flexion distraction with variable electric spring tension allows smooth, counterbalanced flexion, extension or circumduction with the touch of a finger (flexion 17°, extension 14°, side bending 40° total)
- Tilting headpiece with fixed face cushions
- Extending self-locking ankle rest with ankle cuffs
- Upgrade to optional drops with bilateral controls
- Fully modular design for future upgrades
- 23.5" (60 cm) width seamless cushions with multi-density foam
- Lift Capacity: 400 lbs (181 kg)
- Table Length: 69.5" to 76" (177 cm - 193 cm)
- Table Weight: 320 lbs (145 kg)
- Powder coated for durability
- Available in choice of standard upholstery colors

ErgoStyle™ FX Optional Features

- Upgrade to AutoFlexion – In the automatic mode, complete control of manual flexion and reproduction is maintained.
 - All critical information is displayed on the digital interface control:
 - » Range – displays actual degrees of flexion and extension in real time
 - » Speed – variable from 0.5°/second to 6.5°/second
 - » Timer – 1-59 minutes
 - » Display – Flexion degrees, extension degrees, time, speed and angle.
 - Auto Set feature that logs minimum and maximum angles set manually and reproduces in AutoFlexion Mode. Speed, time and angles are adjustable during treatment.
 - Comfort Zone: patient initiated safety feature reduces flexion speed and angle, and gradually increases with each stroke to original set points.
- Patented Accelerator™ III Manual Cocking Drops provide faster acceleration with increased tension. Available for Thoracic, Lumbar and Pelvic sections, and can be added at purchase or as a future upgrade. (Patent #6,679,905)
- Break-away elevating chest with infinite locking capability and single side control assures proper patient "fit"
- Choice of 4 headpiece model upgrades from basic tilting to full drops, all with adjustable face cushions
- T-Bar for manual Flexion
- Upgrade to Designer or Rhino Hyde® Upholstery
- Package includes all air handling hardware, hoses, flow controls and air cylinders with mounting hardware. Does not include air compressor.