

## Eurotech Trademark Table - E9012

### Description

#### Trademark Table Features and Options



---

#### Base Table Standard Features

- Adjustable Face Cushions
- Extension / Flexion Headpiece
- Paper Roller, Cutter and Guide Bar
- Ankle Extension w/ Elevation
- Adjustable Abdominal
- Comfort Cushions

#### Customize Your Table With These Additional Options

- Black Base or White Base
- Elevation (18"-26", 21"-29", or 24"-32")
- Drops: cervical, thoracic, lumbar, pelvic
- Extended Height beyond 26"

#### Standard Table Dimensions

Any table may be ordered from 18" to 26" as measured from the floor to top of cushions.

- Length 65"
- Extended 75"
- Ankle Extension width 24"
- Pelvic cushion width 22"
- Thoracic/Lumbar cushion width 20"
- Cervical cushion width 15"