

ErgoStyle ES2000 Elevation, Adjusting Table - ES9100

Standard Features

- Electric Elevation from 20" to 32" (51cm - 81cm) in 10 seconds with single foot control and expanding covers over lifting mechanism
- Tilting headpiece with fixed face cushions
- 22" (56 cm) width seamless cushions with multi-density foam
- Paper Roll Holder and Cutter
- Extending Ankle Rest
- Powder coated paint
- Available in choice of standard colors

Optional Features

- Upgrade to Auto-Cocking Simple selector switches allow choice of one, two or three drops in any combination for repeated use. With a simple press of the foot control, the table is cocked and ready for use. **Please note: When upgrading to Auto-Cocking, all three drops (thoracic, lumbar and pelvic) must be purchased together. Auto-Cocking on a single drop is not available.*
- Patented modular Accelerator™ II manual and Auto-Cocking drop design provides faster acceleration with increased tension and future upgradeability (Patent # 6,679,905)
- Break-away elevating chest with infinite locking capability and single or dual side control assures proper patient "fit"
- Choice of 4 headpiece model upgrades from basic tilting to full drops, all with adjustable face cushions
- Optional Drops with bilateral controls

Headpiece Options

- ES95000 – Tilt, Fixed Face Cushions
- ES95001 – Tilt, Adjustable Face Cushions
- ES95002 – Tilt, Forward Drop, Adjustable Face Cushions
- ES95003 – Tilt, Elevation, Adjustable Face Cushions
- ES95004 – Tilt, Elevation, Forward Drop, Toggle Drop, Adjustable Face Cushions

Table Options

- ES95005 – Breakaway/Elevating Chest
- ES95006 – Accelerator™ II Thoracic Drop
- ES95007 – Accelerator™ II Lumbar Drop
- ES95008 – Accelerator™ II Pelvic Drop
- ES95009 – Elevating Pelvic Section