



Large Animal Division

by SEDECAL



MIXED ANIMAL PRACTICE SOLUTIONS

SIMPLE. DEPENDABLE. SMART.

SINGLE PANEL DIGITAL IMAGING X-RAY SYSTEM • 14" X 17" FIELD SIZE

800.920.9525
Info@VetRay.com

www.VetRay.com 

Mixed Animal Practice Solutions

Single Panel Digital Imaging X-ray System



Large Animal Division

by SEDECAL

Digital Detector Canon 50G

- Canon Digital flat panel detector with power supply unit
- 14" x 17" full image field size with detachable imaging cable
- Amorphous Silicon Sensor (a-Si)
- GOS scintillator
- 160 micron pixel pitch
- 2.6 million pixels (1464 x 1776)
- 12 bit / 4096 gray scale (14 bits A/D)
- Preview image in 3-5 seconds after exposure, refresh cycle up to 15 seconds
- 7 meter cable on Detector panel to mobile unit



- Canon CXDI software, Version 6.0 for panel operation, image processing & enhancement, and image output
- DICOM 3.0 compatible, including Print Management
- Service Class and Storage Service Class
- Ethernet 10/100 Base T
- Sensor size (W x L x T): 19.3" x 18.8" x 0.9"
- Sensor weight: 5.9 lbs
- Power Supply
- Cables & Multibox

EPX-2

Portable Large Animal Digital X-ray

Canon 50G digital flat panel detector

Lightweight Carrying Case for digital panel, cables and lap top.

Touch Screen
14.1 inch Wide Screen WXGA rugged touchscreen LCD w/Mic

Lap Top Computer

Dell 6420EXR touch screen panel with military grade features for ruggedness and durability

Interface Panel designed for quick setup & takedown



DX-C*

Companion Animal Digital X-ray

EASY AS 1, 2, 3.
Images in 6 seconds.

The Primary Viewing Station can be located anywhere in your network.



Installation and applications training included with purchase.

For more information:

CALL

EMAIL

CLICK



800.920.9525

info@vetray.com

www.vetray.com

*OPTIONS

- SingleTouch Monitor on DX-C
- Portable X-ray Generator
- Mobile Plate Holder
- Mobile Cart